



Cancer Molecular Diagnostics Solid Tumour Service Users' Guide

Version 6
November 2025

Laboratory Contact Details	3
General Enquiries.....	3
Clinical Directors	3
Laboratory	3
Sample selection and dispatch.....	4
Reporting of results	4
Funding	4
Standard of care solid tumour panels.....	5
Lung Mutation Panel.....	5
<i>EGFR</i>	5
<i>ALK</i> Translocation.....	6
<i>ROS1</i> Translocation	6
<i>KRAS</i> Mutation	6
<i>BRAF</i> V600 Mutation.....	7
<i>HER2/Erbb2</i> exon 20 mutation (Research Use Only*)	7
<i>MET</i> exon 14 skipping.....	7
<i>NTRK1/2/3 fusion</i>	7
Lung TKI resistance panel (cfDNA test).....	8
<i>EGFR (T790M and C797S)</i>	8
Colorectal Cancer Panel.....	9
<i>KRAS</i> and <i>NRAS</i> mutation.....	9
<i>BRAF</i> V600 mutation.....	9
<i>ERBB2/HER2</i> copy number gain (Research Use Only*)	9
GIST mutation Panel	9
<i>KIT</i> gene mutation.....	10
<i>PDGFRA</i> gene mutation	10
<i>BRAF</i> V600 mutation.....	10
Melanoma mutation panel.....	11
<i>BRAF</i> V600	11
<i>NRAS</i>	11
<i>KIT</i>	12
Oncomine Focus Library panel (Research Use Only)	12
Gynecological tumors services	
Homologous Recombination Deficiency (HRD) and BRCA testing.....	13
<i>Homologous Recombination Deficiency (HRD) and tumour BRCA testing by Oncomine Comprehensive Assay V3.1</i>	14
<i>BRCA1/2 smMIP (single-molecule molecular inversion probe) panel</i>	15
<i>BRCA1/2 MLPA (Multiplex ligation-dependent probe amplification) large genomic rearrangement detection</i>	15
Appendix 1: Clinically validated mutations for Colorectal, GIST and Melanoma panels.....	16
Appendix 2: Oncomine Focus Hotspot coverage	23
References	34

Laboratory Contact Details

General Enquiries 01 410 3576

Clinical Directors

Prof Stephen Finn 01 416 2998

- Consultant Histopathologist
- Solid Tumour Lead

Dr Niamh Appleby 01 410 3576

- Consultant Haematologist
- Haematology Lead

Laboratory

Dr Cathal O'Brien 01 416 2062

- Chief Clinical Scientist

Dr Georgia Thodi 01 416 6536

- Principal Clinical Scientist

Sarah McCarron 01 410 3576

- Senior Clinical Scientist

Dr Emer Atkinson 01 410 3567

- Senior Clinical Scientist

Sample selection and dispatch

Samples should be submitted as FFPE blocks with an accompanying H&E-stained slide for the following test-requests: standard of care solid tumour tests (Lung Cancer panel, Colorectal Cancer panel, GIST, Melanoma panel) as well as BRCA mutation detection in prostate cancer cases (tBRCA) and HRD status assessment in ovarian cancer cases. FFPE scrolls (2 tubes with 1x10um scroll or 2x10um scrolls in each tube) are requested for the extended Oncomine testing. The H&E-stained slide must be representative of the material in the block so a recent H&E is advised. The slide will be held in the CMD laboratory in case the result needs to be queried in the future, for this reason a slide cut specifically for molecular testing is advised. For germline BRCA testing, a peripheral blood sample should be submitted with the BRCA request form.

Reporting of results

Results are reported by email via SHIRE (CMD LIMS). Reports can be sent to multiple recipients by email, the email addresses for reporting the results should be mentioned in the CMD referral form in the respective sections.

Funding

The National Cancer Control Program currently supports the cost of mutation testing for predictive mutations in lung adenocarcinoma, colorectal adenocarcinoma, malignant melanoma, ovarian cancer, fallopian tube cancer, and primary peritoneal cancer for patients in public hospitals. The cost of testing for private hospitals is not covered under the scope of this funding so testing must be provided under an up-to-date service level agreement.

Tests marked as 'Research Use Only' are not currently covered by NCCP or pharmaceutical company funding and must be paid for by the referring laboratory.

Standard of care Solid tumour service panels

All gene mutations reported in each of the following panels are screened in parallel using the Life Technologies® Oncomine focus assay on the Ion Torrent S5 next generation sequencer.

All samples submitted for routine Lung panel, colorectal panel, melanoma panel and Oncomine testing include NTRK fusion testing by default, as it is considered as tumour agnostic biomarker.

Lung Mutation Panel

The lung mutation panel is focussed on the assessment of treatment guiding *EGFR*, *ERBB2*, *KRAS*, *BRAF* mutations and *ALK*, *ROS1*, *RET*, *NTRK1/2/3* translocations as well as *MET* alterations. In addition to the standard of care assays a number of research use only assays are available singly or as a component of a larger research panel. The panel allows also for the detection of resistance mutations that may emerge upon treatment with TKIs. Particularly *EGFR* (T790M and C797S), *ALK* and *ROS1* hotspot variants inducing resistance to the respective targeted treatments can be identified via lung cancer panel. This information may assist the treating clinician in identifying acquired resistance and selecting the best alternative treatment when necessary.

Sample requirement: FFPE block with accompanying representative H&E-stained slide and histopathology report

Tumour cellularity requirement: >10% tumour cell content (ideally >30%)

Tumour type: Primary or metastatic lesion

Turn around time: 10 working days

EGFR

Mutations of the *EGFR* gene have been associated with an increased likelihood of patient response to therapy targeting the tyrosine kinase domain of *EGFR*. Contrarily, patients with wild-type *EGFR* tend to exhibit worse progression free survival if medicated with *EGFR*-TKIs than a standard chemotherapeutic regimen (Mok et al., 2009). Mutations of the *EGFR* gene vary from point mutations to indels and typically occur in exons 18-21 of the gene. The most frequently observed mutations are associated with a greater likelihood of response to *EGFR*-TKI therapy, however some mutations such as exon 20 insertions may be associated with primary resistance in the erlotinib/gefitinib-naïve setting (Lindeman et al., 2013a; Pirker et al., 2010).

The T790M mutation is thought to confer resistance to *EGFR*-TKI therapy and may be observed in cases of acquired resistance to *EGFR*-TKI therapy (Lindeman et al., 2013a). However, it must also be noted that germline T790M mutations have been observed in a small subset of patients prior to commencement of therapy (Oxnard

et al., 2012). The impact of denovo sensitising and resistance mutations to EGFR-TKIs is currently uncertain, though scant reports suggest that patients will derive benefit from EGFR-TKI therapy (REFERENCE). The link between germline T790M mutation and familial cancer risk is also uncertain and patients with such a mutation profile may be candidates for closer follow up or inclusion in clinical trials.

ALK Translocation

Translocations involving the ALK gene and one of a number of fusion partners are relatively infrequent events in lung cancer and occur in approximately ~5% of cases (Kwak et al., 2010), however, for those patients whose tumours exhibit an ALK translocation significant improvements in progression free survival can be gained by using targeted therapy such as crizotinib (Shaw et al., 2013). The main fusion partner for ALK is the EML4 gene which generates the characteristic EML4-ALK translocation, however, a number of other translocation partners have been identified such as KIF5B (Takeuchi et al., 2009) and TFG (Rikova et al., 2007).

ROS1 Translocation

Translocations of the ROS1 gene are relatively rare in lung cancer and are associated with a response to treatment with the targeted therapy. Fusions of ROS1 protein tyrosine kinase oncogenes with several partner genes have been identified as targetable genetic aberrations in cases of (NSCLC) lacking activating *EGFR*, *KRAS*, *ALK*, *BRAF* or *HER2* oncogene aberrations and it is regarded as the standard of care test for the detection of treatment guiding ROS1 fusions.

RET Translocation

RET gene rearrangements are found in approximately 1–2% of NSCLC cases, predominantly in patients with adenocarcinoma histology, younger age, and minimal or no smoking history. These fusions result in constitutive activation of the RET tyrosine kinase, driving tumor growth and progression. According to the NCCN guidelines, RET fusion-positive NSCLC is best treated with selective RET inhibitors.

KRAS Mutation

Mutations of KRAS are the most commonly observed oncogenic driver mutation in lung adenocarcinoma (Collisson et al., 2014). KRAS mutated tumours with other than KRAS p.(G12C) mutation are not readily treatable with targeted therapies, however, the frequency of KRAS means that it may be a selection criterion for some clinical trials. The presence of a KRAS mutation may, in some cases, be used to infer the absence of another mutation.

BRAF V600 Mutation

Mutations of the BRAF gene occur in a small percentage (<5%) of cases of lung cancer and in contrast to many other cancers which are characterised by BRAF mutations at codon 600, the mutations may be found at other codons. Mutations outside of codon 600 are not currently considered to be clinically relevant. Based on European Medicines Agency approval of combination MEK and BRAF inhibitor therapy for NSCLC we currently report BRAF V600 mutations as part of our standard-of-care lung mutation panel.

HER2/ErbB2 exon 20 mutation (Research Use Only)

Mutations of exon 20 of the *HER2* gene have been observed in approximately 2% of lung adenocarcinomas (Sonobe, Manabe, Wada, & Tanaka, 2006). Although scant, the number of experimental treatment options for patients with such mutations is relatively high. Patients with mutations in exon 20 of the *HER2* gene may be candidates for inclusion in a clinical trial on the basis of their mutation status.

MET exon 14 skipping

Mutations of the *MET* gene which result in alternative splicing with the resultant loss of expression of exon 14 are collectively grouped under the term MET exon 14 skipping. There is growing evidence that METex14 mutation positive NSCLC are generally responsive to MET inhibitor therapy. (Huang, 2020).

NTRK1/2/3 fusion

Fusions of the *NTRK1/2/3* genes occur at a low percentage in many tumours and are also associated with specific cancers such as papillary thyroid cancer, Congenital infantile fibrosarcoma, secretory breast carcinoma and mammary analogue secretory carcinoma (MASC).

The CMD laboratory is currently funded to provide NTRK testing for patients being screened for inclusion in the LOXO-101 trial. In cases where there is a high degree of suspicion that a patient may harbour an NTRK fusion, requests can be forwarded directly to the CMD laboratory for testing.

Lung TKI resistance panel (cfDNA test)

Following treatment with tyrosine kinase inhibitors (TKIs) most patients will suffer relapse due to a number of factors which allow the cancer to develop resistance to targeted therapy. In a significant number of cases the resistance will develop due to an acquired mutation that prevents the TKI from binding to the target molecule. Resistance mutations may be detected in circulating tumour DNA (ctDNA) as a component of the cell free DNA (cfDNA) extracted from plasma samples.

Sample requirement: 2X cell free collection tubes (Roche cfDNA BCT , Streck cfDNA BCT tubes or other equivalent cfDNA collection tubes).

Tumour type: Post-EGFR-TKI therapy*

Turn around time: 10 working days

Sample acceptance: Samples will be accepted for testing Monday to Friday.

Important:

- Always check the **expiry date** on cfDNA tubes before use
- Follow the **manufacturer's instructions** (Roche, Streck, etc) for sample collection and handling
- Store and transport samples at **Room Temperature (18–25°C)**
- Clearly label the **internal and external shipment package** with storage/transport conditions
- **Avoid Friday shipments** – CMD specimen reception operates **Monday to Friday, 08:00–16:00**

*Patients' cancer should have been previously shown to harbour an activating EGFR mutation. Information about the previous mutation should be included with the referral (attached molecular report is advised). A "research only" panel including BRAF and KRAS is available upon request.

EGFR (T790M and C797S)

For patients treated with EGFR-TKI therapy the most common mutation associated with acquired resistance is the EGFR-T790M mutation. Detection of this mutation may be used to select a patient for treatment with a third generation TKI therapy.

For patients with acquired resistance to third generation therapy, resistance may be mediated by MET amplification, HER2 amplification or the development of the EGFR-C797S mutation. The significance of the C797S mutation is thought to differ depending on whether the mutation is encoded in the same gene as the T790M mutation (cis) or on an alternative copy of the gene (trans). Testing using Next Generation Sequencing allows us to report on the cis/trans nature of the mutation which provides valuable information for selecting a patients next line of treatment or selection for clinical trials.

When no resistance mutation is detected in the plasma, confirmation of the negative result in a tissue biopsy is advised.

Colorectal Cancer Panel

The colorectal cancer mutation detection panel is designed to screen for treatment guiding mutations in the *KRAS*, *NRAS* and *BRAF* genes. Regarding *BRAF* gene, only clinically active variants in codon 600 are reported.

Sample requirement: FFPE block with accompanying representative H&E-stained slide and histopathology report

Tumour cellularity requirement: >10% tumour cell content (optimum >30%)

Tumour type: Primary or metastatic lesion

Turn around time: 10 working days

KRAS and *NRAS* mutation

Mutations in codons 12,13,61,117 and 146 of *KRAS* and codons 12,13 and 61 of *NRAS* are known to adversely affect a patient's response to anti-EGFR targeted therapy in colorectal adenocarcinoma (De Roock et al., 2010; Douillard et al., 2013). Only scant data are available for *NRAS* codons 117 and 146 due to the rarity of the mutations and if used to guide therapy a cautious approach would be warranted.

BRAF V600 mutation

Approximately 8% of colorectal cancers will exhibit a *BRAF* V600 mutation. As a general rule *BRAF* V600 mutations tend to be mutually exclusive with mutations of *KRAS* or *NRAS*. *BRAF* mutations tend to be more closely associated with right-sided sessile cancers (Neumann et al., 2012). The combination of *BRAF* mutation status and microsatellite instability status may be used to indicate patient prognosis (Lochhead et al., 2013) or to aid in the diagnosis of Lynch syndrome (Bedeir & Krasinskas, 2011).

The role of *BRAF* as a predictive marker for anti-EGFR therapy is not clearly established and many studies including recent clinical trials have failed to identify a role for *BRAF* mutation in guiding anti-EGFR therapy (Douillard et al., 2013).

GIST mutation Panel

The GIST (GastroIntestinal Tract) cancer mutation detection panel is designed to screen for treatment guiding mutations in the *KIT*, *PDGFRA* and *BRAF* genes. Regarding *BRAF* gene, only clinically active variants in codon 600 are reported.

Sample requirement: FFPE block with accompanying representative H&E-stained slide and histopathology report

Tumour cellularity requirement: >10% tumour cell content (optimum >30%)

Tumour type: Primary or metastatic lesion

Turn around time: 10 working days

***KIT* gene mutation**

Mutations of the *KIT* gene at exons 9,11,13 and 17 are associated with varying degrees of response to a family of tyrosine kinase inhibitors targeting KIT proteins and optimal dose and therapy selection can be exon dependent. Mutations in exon 14 of *KIT* may confer resistance to imatinib and are generally found in the setting of acquired resistance, however, there is some evidence to suggest that lesions with exon 14 mutations may be sensitive to sunitinib (Heinrich et al., 2008).

***PDGFRA* gene mutation**

While less frequent than mutations in *KIT*, mutations in *PDGFRA* are relatively common in GIST, being found in approximately 5% of cases, with mutations tending to be located in exons 12, 14 and 18 of the *PDGFRA* gene (Corless et al., 2005). Mutations in exons 12, 14 and 18 are thought to confer increased sensitivity to therapy targeting mutationally active PDGFRA. Specific exceptions to this rule can be found in the literature, for example, D842 mutations are thought to confer a more resistant phenotype to imatinib and sunitinib (Heinrich et al., 2003, 2008).

***BRAF* V600 mutation**

Mutations in codon 600 of the BRAF gene are a feature of many cancers including GIST however the frequency of mutations is quite low (<1%). While data on the impact of BRAF mutations in GIST is scant, some reports seem to suggest that it may negatively impact the response to imatinib in patients who harbour a concurrent sensitising mutation (Miranda et al., 2012). Secondary BRAF mutations have been suggested as a causal mutation in acquired resistance to imatinib in GIST (Agaram et al., 2008).

GIST responsiveness to agents targeting mutated BRAF protein is currently unknown but patients with a *BRAF* V600 mutation may be eligible for inclusion in a clinical trial on the basis of the result.

Melanoma mutation panel

The Melanoma mutation detection panel is designed to screen for treatment guiding mutations in the *KIT*, *NRAS* and *BRAF* genes. Regarding *BRAF* gene, only clinically active variants in codon 600 are reported.

Melanoma is a highly aggressive form of skin cancer and accounts for a disproportionate number of cancer deaths relative to non-melanoma skin cancer. The mutational landscape of melanoma is dominated by mutations of the *BRAF* gene (Hodis et al., 2012). In addition, *NRAS* mutation is frequently observed and mutations of the *KIT* gene are observed in a number of cases where it may be used to guide therapy.

Sample requirement: FFPE block with accompanying representative H&E-stained slide and histopathology report

Tumour cellularity requirement: >10% tumour cell content (optimum >30%)

Tumour type: Primary or preferably metastatic lesion **

Turn around time: 10 working days

BRAF V600

Mutations of the *BRAF* gene are a common feature of melanoma. The majority of mutations of this gene are found in codon 600 and involve a change at amino acid level from a valine to the glutamic acid (V600E). While the mutation rate of *BRAF* in cutaneous melanoma is thought to be approximately 50%, evidence is starting to accumulate that a lower mutation frequency in the *BRAF* gene is a feature of melanoma in the Irish setting (van den Hurk et al., 2015).

In melanoma patients, response to targeted therapies such as vemurafinib is strongly associated with the presence of a somatic mutation at codon 600 of the *BRAF* gene, therefore the determination of *BRAF* mutation status prior to commencement of therapy is a prescription requirement (da Rocha Dias et al., 2013).

Mutations outside of the V600 codon of *BRAF* are infrequent relative to the dominant V600 mutations and are not used to stratify patients for standard of care therapy.

NRAS

NRAS mutations occur with a frequency of approximately 15% to 26% in melanoma and tend to be mutually exclusive with mutations of the *BRAF* gene (Colombino et al., 2012; Hodis et al., 2012). While *NRAS* mutations are relatively frequent and a likely driver in a subset of melanoma cases, no targeted therapy currently exists for *NRAS*. However, driver mutations in *NRAS* may be independent prognostic factors in melanoma (Jakob et al., 2012) or may be used to select patients for inclusion in clinical trials.

Title: CMD Solid Tumour Service Users Guide

It has also been suggested that NRAS mutation may play a role in therapeutic resistance to imatinib (Hodi et al., 2013) and vemurafinib (Nazarian et al., 2010). Whether the mutations are implicated in acquired or primary resistance to therapy is unclear.

KIT

Mutation and copy number variation of the KIT gene are relatively common occurrences in melanoma at 3% and 6% respectively (Hodis et al., 2012). The melanoma panel has the capacity to detect mutations, but not copy number variations, in a number of coding exons in KIT. Mutations of the KIT gene have been associated with response to treatment with specific tyrosine kinase inhibitors such as imatinib (Carvajal et al., 2011; Hodi et al., 2013).

****In cases of melanoma relapse a biopsy of the recurrent metastatic lesion may be preferable.**

Oncomine Focus Library panel (Research Use Only)

The Oncomine Focus Library Kit assay from Life Technologies is used as the basis for the full research panel. Analysis follows the same clinically validated workflow as for our standard of care assays, however for those assays outside of the scope of our accreditation the same degree of validation cannot be ensured so such assays are designated as 'research use only'.

Panel content (please see note regarding CNV analysis)

DNA part:									
AKT1 (1,3)	EGFR (3,7,12,15,18,19,20,21)	GNA11 (4,5)	JAK3 (11,12,15)	NRAS (2,3,4)					
ALK (21-25)	ERBB2 (8,17-22)	GNAQ (4,5)	KIT (8,9,11,13,17)	PDGFRA (12,14,18)					
AR (6,8)	ERBB3 (2,3,6,8,9)	HRAS (2,3)	KRAS (2,3,4)	PIK3CA (2,5,6,8,10,14,19,21)					
BRAF (11,15)	ERBB4 (18)	IDH1 (4)	MAP2K1 (2,3,6)	RAF1 (7,12)					
CDK4 (2)	ESR1 (9)	IDH2 (4)	MAP2K2 (2)	RET (10,11,13,15,16)					
CTNNB1 (3)	FGFR2 (7,8,9,12,14)	JAK1 (14,15,16)	MET (14,16,19)	ROS1 (36,38)					
DDR2 (5)	FGFR3 (7,9,14,16)	JAK2 (14)	MTOR (30,39,40,43,47,53)	SMO (4,6,8,9)					

DNA part dedicated to CNV:									
ALK	BRAF	CDK4	EGFR	FGFR1	FGFR3	KIT	MET	MYCN	PIK3CA
AR	CCND1	CDK6	ERBB2	FGFR2	FGFR4	KRAS	MYC	PDGFRA	

RNA part:											
ABL1	AKT3	BRAF	ERBB2	ETV1	ETV5	FGFR2	MET	NTRK2	PDGFRA	RAF1	ROS
ALK	AXL	EGFR	ERG	ETV4	FGFR1	FGFR3	NTRK1	NTRK3	PPARG	RET	

Detailed tables provided in appendix 2.

Title: CMD Solid Tumour Service Users Guide

Sample requirement: FFPE scrolls (2 tubes with 1x10um or* 2x10um scrolls in each of them) with accompanying representative H&E-stained slide and histopathology report

Tumour cellularity requirement: >10% tumour cell content (optimum >30%)

Tumour type: Primary or metastatic lesion

Turn around time: 15-20 working days

** depending on tumour volume of the specimen*

Note: Copy number variations (CNV) analysis is not reported as

1. The somatic copy number variant (CNV) analysis requires further standardization and harmonization of the experimental and bio-informatic procedures used. Threshold values for results interpretation have not yet been established.
2. The CNV analysis of the OncoPrint Focus assay (NGS) has not been extensively validated yet.

Retrospective bioinformatic analysis of samples tested per standard of care panel request using the OncoPrint Focus assay is possible for patient selection for clinical trials.

To have samples tested using the OncoPrint gene panel for clinical trials selection please email Prof Stephen Finn (stephen.finn@tcd.ie) Dr Georgia Thodi (GThodi@stjames.ie) or Dr Cathal O'Brien (cpobrien@stjames.ie).

Gynecological tumors services

Homologous Recombination Deficiency (HRD) and BRCA testing.

Homologous recombination deficiency (HRD) is a tumor characteristic that is defined by the inability to accurately repair double-strand breaks (DSBs) in DNA via homologous recombination. HRD can be defined by the presence of germline or somatic variants in BRCA1 or BRCA2, or by the assessment of genomic instability in a tumour sample. BRCA1 and BRCA2 are mutated in ~20% of ovarian cancers, ~10% of HER2 negative breast cancers, and in prostate cancer (~1% BRCA1 and 7-13% BRCA2) (Capoluongo et al., 2017, O'Shaughnessy et al., 2020, Mateo et al., 2024). Additionally, approximately 20% of high-grade ovarian tumours do not have a mutation in either BRCA1 or BRCA2 and are classed as HRD based on genomic instability (Heitz et al 2023).

Patients with a pathogenic/likely pathogenic mutation in BRCA1 or BRCA2 and/or are HRD positive have improved progression free survival when treated with a poly (adenosine diphosphate-ribose) polymerase (PARP) inhibitor.

Pathogenic BRCA1/2 variants are distributed along the entire gene coding regions. Genetic variants include nonsense, frameshift, splice-site, and some missense mutations, as well as large deletions, duplications, and rearrangements.

Current indications for testing for HRD and BRCA mutational status are described on the NCCP website:

<https://www.hse.ie/eng/services/list/5/cancer/profinfo/medonc/sactguidance/testing%20to%20inform%20parp%20inhibitor%20cancer%20treatment%20.html>

For comprehensive HRD status assessment

Sample requirement: FFPE block with accompanying representative H&E-stained slide and histopathology report

Tumour cellularity requirement: >30% tumour*

Tumour type: Pre-chemotherapy biopsy sample (preferred) or post-chemotherapy biopsy sample.

Turn around time: 20-25 working days

* Samples with a lower neoplastic cell content (20-30%) can be analysed but it is highly likely that a GIM (Genomin Instability Metric) will not be reliably generated, thus not reported. BRCA analysis will be reported if the sample meets the minimum internal QC requirements.

Homologous Recombination Deficiency (HRD) and tumour BRCA testing by Oncomine Comprehensive Assay V3.1.

The Oncomine Comprehensive Assay Plus (OCA Plus) enables detection of HRR gene mutations that may cause deficiency in the repair of the repair of DNA double strand breaks via homologous recombination (HRD) as well as reporting the consequence of genomic scarring using the GIM (Genomic Instability Score). The test is intended to identify tumour samples with defects in the homologous recombination repair (HRR) pathway, via the determination of the HRD score [evaluating genomic instability through the assessment of loss of heterozygosity (LOH), telomeric allelic imbalance (TAI), and large-scale state transitions (LST)] and the detection of deleterious variants (nucleotide substitutions, small indels) in the BRCA1 and BRCA2 genes. Copy number variants in the BRCA1 and BRCA2 genes are not reported as CNV analysis has not been fully validated yet as part of the OCA Plus assay.

BRCA germline testing

Peripheral blood (PB) sample requirements: >3ml in EDTA for

BRCA1/2 smMIP (single-molecule molecular inversion probe) panel

Using a smMIP-based in-house designed assay based on probes described by Neveling et al. (Neveling, K. et al. 2017), the open reading frame of BRCA 1 and BRCA2 are sequenced using the Illumina NextSeq550dx platform. The assay is capable of detecting both somatic and germline mutations.

BRCA1/2 MLPA (Multiplex ligation-dependent probe amplification) large genomic rearrangement detection

Large genomic rearrangements are screened for in peripheral blood samples using CE-IVD MLPA probe mixes (MRC-Holland) These copy number changes are often missed by NGS.

Appendix 1: Clinically actionable mutations for Colorectal, GIST and Melanoma panels

Mutations below encompass those mutations that are specifically assessed by the panel. Other, less common mutations may also be detected by the panel but are dependent on panel design. For treatment guiding assays, known SNPs and silent mutations are not reported.

Gene	Exon	Nucleotide Change	Amino acid change
NRAS	4	c.436G>A	p.A146X
NRAS	4	c.436G>T	p.A146X
NRAS	4	c.436G>C	p.A146X
NRAS	3	c.193A>T	p.S65C
NRAS	3	c.183A>T	p.Q61H
NRAS	3	c.183A>T	p.Q61H
NRAS	3	c.183A>G	p.Q61Q
NRAS	3	c.182>183AA>TG	p.Q61L
NRAS	3	c.182>183AA>GG	p.Q61R
NRAS	3	c.181>183CAA>AAG	p.Q61K
NRAS	3	c.182A>G	p.Q61P
NRAS	3	c.182A>T	p.Q61L
NRAS	3	c.182A>G	p.Q61R
NRAS	3	c.181>182CA>TT	p.Q61L
NRAS	3	c.181>182CA>AG	p.Q61R
NRAS	3	c.181C>A	p.Q61K
NRAS	3	c.181C>G	p.Q61E
NRAS	3	c.180>181AC>TA	p.Q61K
NRAS	3	c.179G>A	p.G60E
NRAS	2	c.52G>A	p.A18T
NRAS	2	c.39T>C	p.G13G
NRAS	2	c.39>39GT>TC	p.G13V
NRAS	2	c.38G>A	p.G13D
NRAS	2	c.38G>T	p.G13V
NRAS	2	c.38G>C	p.G13A
NRAS	2	c.37G>C	p.G13R
NRAS	2	c.37G>T	p.G13C
NRAS	2	c.38G>A	p.G13S
NRAS	2	c.36T>C	p.G12G
NRAS	2	c.35G>A	p.G12D

Title: CMD Solid Tumour Service Users Guide

NRAS	2	c.35G>C	p.G12A
NRAS	2	c.35G>T	p.G12v
NRAS	2	c.34>35GG>AA	p.G12N
NRAS	2	c.34G>C	p.G12R
NRAS	2	c.34G>T	p.G12C
NRAS	2	c.34G>A	p.G12S
NRAS	2	c.31G>A	p.A11T
PDGFRA	12	c.1659_1664delGAGGTA	p.R554_Y555delIRY
PDGFRA	12	c.1678_1692del15	p.R560_S564delRVIES
PDGFRA	12	c.1682T>A	p.V561D
PDGFRA	12	c.1694_1695insA	p.S566fs*6
PDGFRA	12	c.1697_1711del15	p.p.S566_E571>K
PDGFRA	12	c.1698_1712del15	p.S566_E571>R
PDGFRA	14	c.1975A>T	p.N659Y
PDGFRA	14	c.1977C>G	p.N659K
PDGFRA	14	c.1977C>A	p.N659K
PDGFRA	18	c.2472C>T	p.V824V
PDGFRA	18	c.2521_2526delAGAGAC	p.R841_D842delRD
PDGFRA	18	c.2524G>T	p.D842Y
PDGFRA	18	c.2524_2525GA>AT	p.D842I
PDGFRA	18	c.2524_2526GAC>TAT	p.D842Y
PDGFRA	18	c.2524_2526delGAC	p.D842delD
PDGFRA	18	c.2524_2532delGACATCATG	p.D842_M844delDIM
PDGFRA	18	c.2524_2535del12	p.I843_D846delIMHD
PDGFRA	18	c.2524_2536>A	p.D842_D846>N
PDGFRA	18	c.2525A>T	p.D842V
PDGFRA	18	c.2526_2537del12	p.I843_D846delIMHD
PDGFRA	18	c.2526_2538>G	p.D842_D846>E
PDGFRA	18	c.2527_2538del12	p.I843_D846delIMHD
PDGFRA	18	c.2528_2539del12	p.I843_S847>T
PDGFRA	18	c.2530_2541del12	p.M844_S847delMHDS
PDGFRA	18	c.2536G>T	p.D846Y
KIT	9	c.1509_1510insGCCTAT	p.Y503_F504insAY
KIT	9	c.1526A>T	p.K509I
KIT	11	c.1651_1665del15	p.P551_V555del
KIT	11	c.1652C>T	p.P551L
KIT	11	c.1653_1670del18	p.M552_W557del
KIT	11	c.1654A>C	p.M552L
KIT	11	c.1654_1659delATGTAT	p.M552_Y553del

Title: CMD Solid Tumour Service Users Guide

KIT	11	c.1654_1662del9	p.M552_E554del
KIT	11	c.1654_1668del15	p.M552_Q556del
KIT	11	c.1654_1671del18	p.M552_W557del
KIT	11	c.1655_1660delTGATG	p.M552_E554>K
KIT	11	c.1655_1666del12	p.M552_Q556>K
KIT	11	c.1655_1672del18	p.M552_W557del
KIT	11	c.1656_1670del15	p.Y553_W557del
KIT	11	c.1656_1673del18	p.Y553_K558
KIT	11	c.1657T>A	p.Y553N
KIT	11	c.1657_1668del12	p.Y553_Q556del
KIT	11	c.1657_1674del18	p.Y553_K558del
KIT	11	c.1658_1720del63	p.Y553_T574>S
KIT	11	c.1660G>A	p.E554K
KIT	11	c.1660_1674del15	p.E554_K558del
KIT	11	c.1661_1675del15	p.E554_K558del
KIT	11	c.1662_1685del24	p.V555_E562del
KIT	11	c.1663G>A	p.V555I
KIT	11	c.1663_1668delGTACAG	p.V555_Q556del
KIT	11	c.1663_1674del12	p.V555_K558del
KIT	11	c.1663_1677del15	p.V555_V559del
KIT	11	c.1663_1680del18	p.V555_V560del
KIT	11	c.1663_1713del51	p.V555_I571del
KIT	11	c.1663_1719del57	p.V555_P573del
KIT	11	c.??>?del?	p.Q556>K558del
KIT	11	c.1666_1680del15	p.Q556_V560del
KIT	11	c.1666_1728del63	p.Q556_L576del
KIT	11	c.1667_1672delAGTGGA	p.W557_K558del
KIT	11	c.1667_1681del15	p.W557_E561del
KIT	11	c.1668_1673delGTGGAA	p.W557_K558del
KIT	11	c.1668_1679del12	p.Q556_V560>H
KIT	11	c.1669T>A	p.W557R
KIT	11	c.1669T>C	p.W557R
KIT	11	c.1669T>G	p.W557G
KIT	11	c.1669_1671delTGG	p.W557del
KIT	11	c.1669_1672TGGGA>G	p.W557_K558>E
KIT	11	c.1669_1674delTGGGAAG	p.W557_K558del
KIT	11	c.1669_1674delTGGGAAG	p.W557_K558del
KIT	11	c.1669_1677del9	p.W557_V559del
KIT	11	c.1669_1680del12	p.W557_V560del

Title: CMD Solid Tumour Service Users Guide

KIT	11	c.1669_1683del15	p.W557_E561del
KIT	11	c.1670G>C	p.W557S
KIT	11	c.1670_1675delGGAAGG	p.W557_V559>F
KIT	11	c.1670_1717del48	p.W557
KIT	11	c.1671_1673GAA>TCC	p.W557_K558>CP
KIT	11	c.1671_1676delGAAGGT	p.W557_V559>C
KIT	11	c.1671_1679del9	p.W557_V560>C
KIT	11	c.1672A>G	p.K558E
KIT	11	c.1672_1676AAGGT>TCTTC	p.K558_V559>SS
KIT	11	c.1672_1677delAAGGTT	p.K558_V559del
KIT	11	c.1672_1680del9	p.K558_V560del
KIT	11	c.1672_1686del15	p.K558_E562del
KIT	11	c.1672_1692del21	p.K558_N564del
KIT	11	c.1673A>G	p.K558R
KIT	11	c.1673_1678delAGGTTG	p.K558>V560>I
KIT	11	c.1673_1687del15	p.K558>E562del
KIT	11	c.1673_1674insTCC	p.K558>NP
KIT	11	c.1674G>A	p.K558K
KIT	11	c.1674_1674G>TCCT	p.K558>NP
KIT	11	c.1674_1679delGGTTGT	p.K558>V560>N
KIT	11	c.1675G>A	p.V559I
KIT	11	c.1675_1677delGTT	p.V559del
KIT	11	c.1675_1680delGTTGTT	p.V559_V560del
KIT	11	c.1675_1683del9	p.V559_E561del
KIT	11	c.1675_1695del21	p.V559_G565del
KIT	11	c.1676T>A	p.V559D
KIT	11	c.1676T>G	p.V559G
KIT	11	c.1676T>C	p.V559A
KIT	11	c.1676_1684del9	p.V559_E561del
KIT	11	c.1678_1680delGTT	p.V560del
KIT	11	c.1678_1728del51	p.V560_L576del
KIT	11	c.1679T>A	p.V560D
KIT	11	c.1679T>G	p.V560G
KIT	11	c.1679_1680TT>AG	p.V560E
KIT	11	c.1679_1681delTTG	p.V560del
KIT	11	c.1681G>A	p.E561K
KIT	11	c.1683G>A	p.E561E
KIT	11	c.1684G>A	p.E562K
KIT	11	c.1687_1728del42	p.I563_L576del

Title: CMD Solid Tumour Service Users Guide

KIT	11	c.1690_1728del39	p.N564_L576del
KIT	11	c.1690_1734del45	p.N564_Y578del
KIT	11	c.1696A>G	p.N566D
KIT	11	c.1698C>T	p.N566N
KIT	11	c.1701T>A	p.N567K
KIT	11	c.1702T>G	p.Y568D
KIT	11	c.1702_1722del21	p.Y568_T574del
KIT	11	c.1708_1728del21	p.Y570_L576del
KIT	11	c.1711_1728del18	p.I571_L576del
KIT	11	c.1714G>A	p.D572N
KIT	11	c.1718C>T	p.P573L
KIT	11	c.1726C>T	p.L576F
KIT	11	c.1726>1728delCTT	p.L576del
KIT	11	c.1727T>C	p.L576P
KIT	11	c.1729C>T	p.P577S
KIT	11	c.1735_1737delGAT	p.D579del
KIT	11	c.1745G>A	p.W582*
KIT	11	c.1746G>A	p.W582*
KIT	11	c.1751T>C	p.F584S
KIT	11	c.1755C>T	p.P585P
KIT	13	c.1924A>G	p.K642E
KIT	13	c.1928T>C	p.V643A
KIT	13	c.1961T>C	p.V654A
KIT	14	c.2008>2009AC>GA	p.T670E
KIT	14	c.2009C>T	p.T670I
KIT	17	c.2433T>C	p.F811F
KIT	17	c.2446G>T	p.D816Y
KIT	17	c.2446G>C	p.D816H
KIT	17	c.2446_2447GA>TT	p.D816F
KIT	17	c.2446_2447GA>AT	p.D816I
KIT	17	c.2447A>G	p.D816G
KIT	17	c.2447A>T	p.D816V
KIT	17	c.2448C>G	p.D816E
KIT	17	c.2453A>G	p.K818R
KIT	17	c.2458G>T	p.D820Y
KIT	17	c.2458G>C	p.D820H
KIT	17	c.2459A>G	p.D820G
KIT	17	c.2459A>T	p.D820V
KIT	17	c.2460T>A	p.D820E

Title: CMD Solid Tumour Service Users Guide

KIT	17	c.2464A>T	p.N822Y
KIT	17	c.2466T>A	p.N822K
KIT	17	c.2466T>G	p.N822K
KIT	17	c.2467T>G	p.Y823D
KIT	17	c.2468A>G	p.Y823C
KIT	17	c.2473G>A	p.V825I
KIT	17	c.2474T>C	p.V825A
BRAF	15	c.1799_1801delTGA	p.V600_K601>E
BRAF	15	c.1799_1800TG>AA	p.V600E
BRAF	15	c.1799_1800TG>AT	p.V600D
BRAF	15	c.1799T>C	p.V600A
BRAF	15	c.1799T>A	p.V600E
BRAF	15	c.1799T>G	p.V600G
BRAF	15	c.1798_1977GT>CA	p.V600Q
BRAF	15	c.1798_1799GT>AA	p.V600K
BRAF	15	c.1798_1799GT>AG	p.V600R
BRAF	15	c.1797_1799AGT>GAG	p.V600R
BRAF	15	c.1798G>A	p.V600M
BRAF	15	c.1798G>C	p.V600L
BRAF	15	c.1798G>T	p.V600L
KRAS	4	c.437C>T	p.A146V
KRAS	4	c.436G>A	p.A146T
KRAS	4	c.436G>C	p.A146P
KRAS	4	c.351A>C	p.K117N
KRAS	4	c.351A>T	p.K117N
KRAS	3	c.183A>C	p.Q61H
KRAS	3	c.183A>T	p.Q61H
KRAS	3	c.182A>C	p.Q61P
KRAS	3	c.182A>G	p.Q61R
KRAS	3	c.182A>T	p.Q61L
KRAS	3	c.181C>A	p.Q61K
KRAS	3	c.181C>G	p.Q61E
KRAS	3	c.180_181TC>AA	p.Q61K
KRAS	3	c.176C>G	p.A59G
KRAS	3	c.176C>A	p.A59E
KRAS	3	c.175G>A	p.A59T
KRAS	2	c.39C>G	p.G13G
KRAS	2	c.39C>T	p.G13G
KRAS	2	c.39C>A	p.G13G

Title: CMD Solid Tumour Service Users Guide

KRAS	2	c.38_39GG>TT	p.G13V
KRAS	2	c.38_39GC>TG	p.G13V
KRAS	2	c.38_39GC>AT	p.G13D
KRAS	2	c.38_39GC>AA	p.G13E
KRAS	2	c.37_39GGC>CGT	p.G13R
KRAS	2	c.38G>A	p.G13D
KRAS	2	c.38G>C	p.G13A
KRAS	2	c.38G>T	p.G13V
KRAS	2	c.37G>T	p.G13C
KRAS	2	c.37G>A	p.G13S
KRAS	2	c.37G>C	p.G13R
KRAS	2	c.36_37TG>AT	p.G13C
KRAS	2	c.36_37insGGT	p.G12>G13insG
KRAS	2	c.36T>C	p.G12G
KRAS	2	c.36T>A	p.G12G
KRAS	2	c.35_36GT>AC	p.G12D
KRAS	2	c.36_36GT>TC	p.G12V
KRAS	2	c.35_36GT>AA	p.G12E
KRAS	2	c.34_36GGT>TGG	p.G12W
KRAS	2	c.34_36GGT>TGC	p.G12C
KRAS	2	c.35G>T	p.G12V
KRAS	2	c.35G>A	p.G12D
KRAS	2	c.35G>C	p.G12A
KRAS	2	c.34_35GG>TA	p.G12Y
KRAS	2	c.34_35GG>AT	p.G12I
KRAS	2	c.34_35GG>TT	p.G12F
KRAS	2	c.34_35GG>CT	p.G12L
KRAS	2	c.34G>T	p.G12C
KRAS	2	c.34G>A	p.G12S
KRAS	2	c.34G>C	p.G12R

Appendix 2: Oncomine Focus Hotspot coverage

Mutation hotspots

Gene	Exon	Amino Acids	Accession	Chromosome
AKT1	3	17-51	NM_005163	14
ALK	11	660-680	NM_004304	2
ALK	21	1121-1150	NM_004304	2
ALK	22	1151-1171	NM_004304	2
ALK	23	1173-1214	NM_004304	2
ALK	24	1217-1247	NM_004304	2
ALK	25	1252-1277	NM_004304	2
ALK	27	1334-1357	NM_004304	2
ALK	29	1584-1620	NM_004304	2
ALK	29	1499-1541	NM_004304	2
ALK	29	1428-1470	NM_004304	2
BRAF	5	204-228	NM_004333	7
BRAF	6	249-273	NM_004333	7
BRAF	7	288-309	NM_004333	7
BRAF	8	328-349	NM_004333	7
BRAF	10	418-438	NM_004333	7
BRAF	11	451-477	NM_004333	7
BRAF	13	594-565	NM_004333	7
BRAF	13	507-513	NM_004333	7
BRAF	15	582-610	NM_004333	7
BRAF	18	736-766	NM_004333	7
CDK4	2	8-49	NM_000075	12
CDK4	5	190-210	NM_000075	12
CDK4	6	212-224	NM_000075	12
CDK4	7	268-273	NM_000075	12
CDK4	8	286-302	NM_000075	12
CTNNB1	3	10-47	NM_001904	3
DDR2	5	93-134	NM_006182	1
EGFR	3	97-123	NM_005228	7
EGFR	7	280-296	NM_005228	7
EGFR	12	474-499	NM_005228	7
EGFR	15	575-601	NM_005228	7
EGFR	18	695-725	NM_005228	7
EGFR	19	729-761	NM_005228	7
EGFR	20	762-799	NM_005228	7
EGFR	21	827-865	NM_005228	7

Title: CMD Solid Tumour Service Users Guide

ERBB2	8	302-326	NM_004448	17
ERBB2	17	662-690	NM_004448	17
ERBB2	18	704-736	NM_004448	17
ERBB2	19	753-769	NM_004448	17
ERBB2	20	770-796	NM_004448	17
ERBB2	21	840-881	NM_004448	17
ERBB2	22	884-904	NM_004448	17
ERBB2	24	969-989	NM_004448	17
ERBB2	25	1010-1050	NM_004448	17
ERBB3	2	47-77	NM_001982	12
ERBB3	3	88-120	NM_001982	12
ERBB3	6	206-244	NM_001982	12
ERBB3	8	293-312	NM_001982	12
ERBB3	9	331-369	NM_001982	12
ERBB4	18	694-713	NM_005235	2
ESR1	8	533-572	NM_00125	6
FGFR2	7	251-274	NM_000141	10
FGFR2	8	315-348	NM_000141	10
FGFR2	9	363-398	NM_000141	10
FGFR2	12	546-557	NM_000141	10
FGFR2	14	631-662	NM_000141	10
FGFR3	3	80-121	NM_000142	4
FGFR3	7	248-277	NM_000142	4
FGFR3	9	368-402	NM_000142	4
FGFR3	14	632-653	NM_000142	4
FGFR3	16	691-719	NM_000142	4
GNA11	4	168-201	NM_002067	19
GNA11	5	203-219	NM_002067	19
GNAQ	4	166-195	NM_002072	9
GNAQ	5	206-245	NM_002072	9
HRAS	2	1-25	NM_005343	11
HRAS	3	43-81	NM_005343	11
IDH1	4	101-134	NM_005896	2
IDH2	4	134-176	NM_002168	15
JAK1	14	638-662	NM_002227	1
JAK1	15	686-705	NM_002227	1
JAK1	16	706-742	NM_002227	1
JAK2	14	604-621	NM_004972	9
JAK3	11	484-507	NM_000215	19
JAK3	12	525-567	NM_000215	19
JAK3	15	641-681	NM_000215	19

Title: CMD Solid Tumour Service Users Guide

KIT	8	412-433	NM_000222	4
KIT	9	495-513	NM_000222	4
KIT	10	526-549	NM_000222	4
KIT	11	550-587	NM_000222	4
KIT	13	628-660	NM_000222	4
KIT	17	792-828	NM_000222	4
KRAS	2	6-36	NM_004985	12
KRAS	3	38-65	NM_004985	12
KRAS	4	114-149	NM_004985	12
MAP2K1	2	44-81	NM_002755	15
MAP2K1	3	106-144	NM_002755	15
MAP2K1	6	191-225	NM_002755	15
MAP2K2	2	32-67	NM_030662	19
MET	2	66-102	NM_000245	7
MET	2	360-391	NM_000245	7
MET	11	799-837	NM_000245	7
MET	14	982-1009	NM_000245	7
MET	15	1063-1086	NM_000245	7
MET	16	1088-1113	NM_000245	7
MET	19	1229-1255	NM_000245	7
MTOR	30	1457-1489	NM_004958	1
MTOR	39	1789-1810	NM_004958	1
MTOR	40	1872-1904	NM_004958	1
MTOR	43	1999-2011	NM_004958	1
MTOR	43	1971-1979	NM_004958	1
MTOR	47	2187-2221	NM_004958	1
MTOR	53	2394-2434	NM_004958	1
NRAS	2	4-30	NM_002524	1
NRAS	3	43-68	NM_002524	1
NRAS	4	113-149	NM_002524	1
PDGFRA	12	552-583	NM_006206	4
PDGFRA	14	645-667	NM_006206	4
PDGFRA	18	820-853	NM_006206	4
PIK3CA	2	24-60	NM_006218	3
PIK3CA	2	79-113	NM_006218	3
PIK3CA	5	317-351	NM_006218	3
PIK3CA	6	354-377	NM_006218	3
PIK3CA	8	418-433	NM_006218	3
PIK3CA	8	450-468	NM_006218	3
PIK3CA	10	533-554	NM_006218	3

Title: CMD Solid Tumour Service Users Guide

PIK3CA	14	692-729	NM_006218	3
PIK3CA	19	899-923	NM_006218	3
PIK3CA	21	1018-1050	NM_006218	3
RAF1	7	246-273	NM_002880	3
RAF1	12	407-441	NM_002880	3
RET	10	609-626	NM_020630	10
RET	11	628-653	NM_020630	10
RET	13	763-785	NM_020630	10
RET	15	876-909	NM_020630	10
RET	16	911-923	NM_020630	10
ROS1	36	1942-1980	NM_002944	6
ROS1	38	2003-2036	NM_002944	6
SMO	6	392-418	NM_005631	7
SMO	8	454-488	NM_005631	7
SMO	9	512-542	NM_005631	7

Fusion Genes

Fusion Gene	Fusion	Exons	Cosmic Fusion number
ABL1	EML1-ABL1	E17A2	
AKT3	MAGI3-AKT3	M10A3	
ALK	ACTG2-ALK	A2A18	
ALK	EML4-ALK	E13A20	AB462411
ALK	EML4-ALK	E18A20	COSF487
ALK	EML4-ALK	E20A20	COSF409
ALK	EML4-ALK	E2A20	COSF478
ALK	EML4-ALK	E6A17	
ALK	EML4-ALK	E6A18	
ALK	EML4-ALK	E6A19	COSF1296
ALK	EML4-ALK	E6aA20	AB374361
ALK	GTF2IRD1-ALK	G7A20	
ALK	HIP1-ALK	H21A20	
ALK	HIP1-ALK	H28A20	
ALK	KIF5B-ALK	K15A20	COSF1381
ALK	KIF5B-ALK	K17A20	COSF1257
ALK	KIF5B-ALK	K24A20	COSF1058
ALK	KLC1-ALK	K9A20	COSF1276
ALK	MEMO1-ALK	M2A7	
ALK	STRN-ALK	S3A20	COSF1430
ALK	TFG-ALK	T4A20	COSF424
ALK	TFG-ALK	T5A20	COSF426

Title: CMD Solid Tumour Service Users Guide

ALK	TFG-ALK	T6A20	COSF428
ALK	TPM1-ALK	T8A20	
ALK	TPM3-ALK	T7A20	
ALK	TPR-ALK	T15A20	
ALK	VCL-ALK	V16A20	COSF1057
ALK	A2M-ALK	A22A19	
ALK	ATIC-ALK	A7A20	
ALK	CARS-ALK	C17A20	
ALK	CLTC-ALK	C31A20	
ALK	DCTN1-ALK	D26A20	
ALK	EML4-ALK	E21A20	
ALK	RANBP2-ALK	R18A20	
ALK	SEC31L1_SEC31A-ALK	S21A20	
ALK	SEC31L1_SEC31A-ALK	S22A20	
ALK	TRAF1-ALK	T6A20	
ALK	TPM4-ALK	T7A20	
ALK	C2orf44-ALK	C4A20	
ALK	EML4-ALK	E13A20	COSF1062
ALK	EML4-ALK	E14A20	COSF1064
ALK	EML4-ALK	E14A20	COSF477
ALK	EML4-ALK	E17A20	COSF1366
ALK	EML4-ALK	E17A20	COSF1367
ALK	EML4-ALK	E20A20	COSF730
ALK	EML4-ALK	E2A20	COSF479
ALK	EML4-ALK	E6bA20	AB374362
ALK	EML4-ALK	E7A20	
ALK	KIF5B-ALK	K15A20	COSF1060
ALK	NCOA1-ALK	N21A1	
ALK	EML4-ALK	E3p53insA20	
ALK	SMEK2-ALK	S9A2	
ALK	EML4-ALK	E15A20	COSF413
ALK	PRKAR1A-ALK	P2A20	
ALK	EML4-ALK	E6ins18A20	
ALK	EML4-ALK	E13ins90A20	
ALK	EML4-ALK	E14ins2del52A20	
ALK	EML4-ALK	E14ins124A20	
ALK	EML4-ALK	E14del36A20	
ALK	EML4-ALK	E17ins65A20	
ALK	EML4-ALK	E17ins68A20	
ALK	EML4-ALK	E17ins30A20_V8a	
ALK	CLIP4-ALK	C12A23	

Title: CMD Solid Tumour Service Users Guide

BRAF	AKAP9-BRAF	A8B9	COSF1013.1
BRAF	CDC27-BRAF	C16B9	
BRAF	FAM131B-BRAF	F2B9	COSF1189
BRAF	FCHSD1-BRAF	F13B9	COSF404
BRAF	KIAA1549-BRAF	K16B10	COSF509
BRAF	KIAA1549-BRAF	K17B10	
BRAF	KIAA1549-BRAF	K18B9	
BRAF	PAPSS1-BRAF	P5B9	
BRAF	SLC45A3-BRAF	S1B8	COSF871
BRAF	SND1-BRAF	S16B9	
BRAF	TRIM24-BRAF	T9B9	
BRAF	KIAA1549-BRAF	K15B9	
BRAF	KIAA1549-BRAF	K15B11	
BRAF	TAX1BP1-BRAF	T8B11	
BRAF	AGTRAP-BRAF	A5B8	COSF828
EGFR	EGFR-EGFR	E1E8	
ERBB2	WIPF2-ERBB2	W1E4	
ERG	SLC45A3-ERG	S1E4	COSF1138
ERG	TMPRSS2-ERG	T1E2	COSF23
ERG	TMPRSS2-ERG	T1E3	COSF24
ERG	TMPRSS2-ERG	T1E4	COSF38
ERG	TMPRSS2-ERG	T1E5	COSF26
ERG	TMPRSS2-ERG	T1E6	COSF36
ERG	TMPRSS2-ERG	T1EIIIc_4	
ERG	TMPRSS2-ERG	T2E2	COSF27
ERG	TMPRSS2-ERG	T2E4	COSF28
ERG	TMPRSS2-ERG	T2EIIIc_4	
ERG	TMPRSS2-ERG	T2E5	COSF29
ERG	TMPRSS2-ERG	T3E4	COSF30
ERG	TMPRSS2-ERG	T4E4	COSF18
ERG	TMPRSS2-ERG	T4E5	COSF17
ERG	TMPRSS2-ERG	T5E4	COSF16
ERG	TMPRSS2-ERG	T5E5	
ETV1	TMPRSS2-ETV1	T1E5	
ETV1	TMPRSS2-ETV1	T2E5	
ETV1	TMPRSS2-ETV1	T1E4	
ETV1	TMPRSS2-ETV1	T1bE4	
ETV4	TMPRSS2-ETV4	T1bE3	
ETV5	TMPRSS2-ETV5	T1E2	EU314929
ETV5	TMPRSS2-ETV5	T3E2	EU314930
ETV5	TMPRSS2-ETV5	T1bE2	

Title: CMD Solid Tumour Service Users Guide

FGFR1	BAG4-FGFR1	B2F6	
FGFR1	ERLIN2-FGFR1	E8F2	
FGFR1	BAG4-FGFR1	B1F8	
FGFR1	FGFR1-TACC1	F17T7	
FGFR2	FGFR2-AFF3	F17A8	
FGFR2	FGFR2-BICC1	F17B2	
FGFR2	FGFR2-CASP7	F17C2	
FGFR2	FGFR2-CIT	F17C23	
FGFR2	SLC45A3-FGFR2	S1F1	
FGFR2	SLC45A3-FGFR2	S1F2	
FGFR2	FGFR2-KIAA1967_CCAR2	F17C4	
FGFR2	FGFR2-MGEA5	F17M12	
FGFR2	FGFR2-OFD1	F17O3	
FGFR2	FGFR2-TACC3	F17T11	
FGFR3	FGFR3-AES	F17A2	
FGFR3	FGFR3-BAIAP2L1	F17B2	COSF1347
FGFR3	FGFR3-ELAVL3	F17E2	
FGFR3	FGFR3-TACC3	F15T11	
FGFR3	FGFR3-TACC3	F16T10	COSF1359
FGFR3	FGFR3-TACC3	F16T11	COSF1348
FGFR3	FGFR3-TACC3	F17T10	COSF1434
FGFR3	FGFR3-TACC3	F17T11	
FGFR3	FGFR3-TACC3	F17T5	
FGFR3	FGFR3-TACC3	F17T6	
FGFR3	FGFR3-TACC3	F17T8	
FGFR3	FGFR3-TACC3	F17T9	
FGFR3	FGFR3-TACC3	F17T13	
FGFR3	FGFR3-TACC3	F18T7	
FGFR3	FGFR3-TACC3	F14T11	
FGFR3	FGFR3-TACC3	F18T4and5	
FGFR3	FGFR3-TACC3	F18T10	
FGFR3	FGFR3-TACC3	F18T10	
FGFR3	FGFR3-TACC3	F18T11	
FGFR3	FGFR3-TACC3	TruncatedF17T4	
FGFR3	FGFR3-TACC3	F17T7	
FGFR3	FGFR3-TACC3	F17T10	
FGFR3	FGFR3-TACC3	F17T11	
FGFR3	FGFR3-TACC3	F17T11	
FGFR3	FGFR3-TACC3	F17intron17T4	
FGFR3	FGFR3-TACC3	F17Intron17T9	
MBIP	AXL-MBIP	A20M4	

Title: CMD Solid Tumour Service Users Guide

MET	BAIAP2L1-MET	B9M15	
MET	C8orf34-MET	C2M15	
MET	CAPZA2-MET	C4M11	
MET	MET-MET	M13M15	
MET	OXR1-MET	O9M13	
MET	PTPRZ1-MET	P1M2	
MET	PTPRZ1-MET	P3M2	
MET	PTPRZ1-MET	P8M2	
MET	TFG-MET	T5M15	
MET	TPR-MET	T4M15	
NTRK1	NTRK1-DYNC2H1	N17D85	
NTRK1	MPRIP-NTRK1	M14N12	
NTRK1	MPRIP-NTRK1	M18N12	
NTRK1	MPRIP-NTRK1	M21N12	
NTRK1	SQSTM1-NTRK1	S5N10	
NTRK1	SSBP2-NTRK1	S12N12	
NTRK1	TFG-NTRK1	T6N10	
NTRK1	TPM3-NTRK1	T7N10	COSF1318
NTRK1	TPR-NTRK1	T21N10	
NTRK1	BCAN-NTRK1	B13N11	
NTRK1	LMNA-NTRK1	L2N11	
NTRK1	NFASC-NTRK1	N20N10	
NTRK1	RNF213-NTRK1	R15N12	
NTRK1	CD74-NTRK1	C7N10	
NTRK1	CEL-NTRK1	C7N7	
NTRK1	IRF2BP2-NTRK1	I1N10	
NTRK1	TPR-NTRK1	T6N12	
NTRK1	TPR-NTRK1	T6N12	
NTRK1	TPR-NTRK1	T21N9	
NTRK2	AFAP1-NTRK2	A14N12	
NTRK2	AGBL4-NTRK2	A6N16	
NTRK2	NACC2-NTRK2	N4N13	
NTRK2	QKI-NTRK2	Q6N16	
NTRK2	SQSTM1-NTRK2	S5N17	
NTRK2	TRIM24-NTRK2	T12N15	
NTRK2	VCL-NTRK2	V16N12	
NTRK3	ETV6-NTRK3	E4N15	COSF823
NTRK3	ETV6-NTRK3	E5N15	COSF571
NTRK3	BTBD1-NTRK3	B4N14	
NTRK3	COX5A-NTRK3	C1N15	
NTRK3	ETV6-NTRK3	E4N14	

Title: CMD Solid Tumour Service Users Guide

NTRK3	ETV6-NTRK3	E5N14	
PDGFRA	SCAF11-PDGFRA	S1P2	
PPARG	PAX8-PPARG	P10P2	COSF1219
PPARG	PAX8-PPARG	P7P2	COSF1223
PPARG	PAX8-PPARG	P8P2	COSF1215
PPARG	PAX8-PPARG	P9P2	COSF1217
PTPN3	ALK-PTPN3	A11P3	
RAF1	B4GALT1-RAF1	B1R8	
RAF1	ESRP1-RAF1	E13R6	COSF826
RET	CCDC6-RET	C1R12	COSF1271
RET	CUX1-RET	C10R12	
RET	ERC1_ELKS-RET	E11R12	COSF1508
RET	ERC1-RET	E12R12	
RET	ERC1-RET	E17R12	
RET	ERC1-RET	E7R12	
RET	FKBP15-RET	F25R12	
RET	GOLGA5-RET	G7R12	
RET	HOOK3-RET	H11R12	
RET	KIAA1468-RET	K10R12	
RET	KIF5B-RET	K15R12	COSF1232
RET	KIF5B-RET	K16R12	COSF1230
RET	KIF5B-RET	K22R12	COSF1253
RET	KIF5B-RET	K23R12	COSF1234
RET	KIF5B-RET	K24R11	COSF1262
RET	KIF5B-RET	K24R8	COSF1236
RET	NCOA4-RET	N7R12	
RET	PCM1-RET	P29R12	
RET	PRKAR1A-RET	P7R12	
RET	TBL1XR1-RET	T9R12	
RET	TRIM27-RET	T3R12	
RET	TRIM33-RET	T16R12	
RET	ACBD5-RET	A11R12	
RET	AFAP1-RET	A3R12	
RET	AKAP13-RET	A35R12	
RET	NCOA4_ELE1-RET	E6R12	
RET	SPECC1L-RET	S10R12	
RET	KTN1-RET	K29R12	COSF1513
RET	RUFY2-RET	R9R12	
RET	KIF5B-RET	K15R11	COSF1255
RET	SPECC1L-RET	S10R11	
RET	TBL1XR1-RET	T9R11	

Title: CMD Solid Tumour Service Users Guide

RET	TRIM24-RET	T9R12	
RET	AKAP13-RET	A36R12	
RET	CCDC6-RET	C1R13	
RET	CCDC6-RET	C2R12	
RET	CCDC6-RET	C8R11	
RET	CCDC6-RET	C8R12	
RET	CCDC6-RET	C1R12	
RET	CCDC6-RET	C1R11	
RET	CCDC6-RET	C1R11	
ROS1	CD74-ROS1	C6R32	COSF1202
ROS1	CD74-ROS1	C6R34	COSF1200
ROS1	CEP85L-ROS1	C8R36	
ROS1	EZR-ROS1	E10R34	COSF1267
ROS1	GOPC-ROS1	G4R36	COSF1188
ROS1	GOPC-ROS1	G8R35	COSF1139
ROS1	LRIG3-ROS1	L16R35	COSF1269
ROS1	SDC4-ROS1	S2R32	COSF1265
ROS1	SDC4-ROS1	S2R34	
ROS1	SDC4-ROS1	S4R32	COSF1278
ROS1	SDC4-ROS1	S4R34	COSF1280
ROS1	SLC34A2-ROS1	S4R32	COSF1197
ROS1	SLC34A2-ROS1	S4R34	COSF1198
ROS1	CLTC-ROS1	C31R35	
ROS1	ERC1-ROS1	E11R36	
ROS1	HLA_A-ROS1	H7R34	
ROS1	KIAA1598-ROS1	K11R36	
ROS1	MYO5A-ROS1	M23R35	
ROS1	PPFIBP1-ROS1	P9R35	
ROS1	PWWP2A-ROS1	P1R36	
ROS1	TPM3-ROS1	T3R36	
ROS1	ZCCHC8-ROS1	Z2R36	
ROS1	TFG-ROS1	T4R35	
ROS1	CD74-ROS1	C6R35	
ROS1	EZR-ROS1	E10R35	
ROS1	MSN-ROS1	M9R34	
ROS1	CCDC6-ROS1	C5R35	
ROS1	CD74-ROS1	C4R33	
ROS1	KDEL2-ROS1	K5R35	
ROS1	SLC34A2-ROS1	S13R32	COSF1259
ROS1	SLC34A2-ROS1	S13R34	COSF1261
ROS1	SLC34A2-ROS1	S13R36	

Title: CMD Solid Tumour Service Users Guide

ROS1	TPM3-ROS1	T7R35	COSF1273
ROS1	CLIP1-ROS1	C19R36	
ROS1	KDEL2-ROS1	K5Rintron34	

References

- Agaram, N. P., Wong, G. C., Guo, T., Maki, R. G., Singer, S., DeMatteo, R. P., ... Antonescu, C. R. (2008). Novel V600E BRAF mutations in imatinib-naïve and imatinib-resistant gastrointestinal stromal tumors. *Genes Chromosomes and Cancer*, 47(10), 853–859. doi:10.1002/gcc.20589
- Bedeir, A., & Krasinskas, A. M. (2011). Molecular diagnostics of colorectal cancer. *Archives of Pathology & Laboratory Medicine*, 135(5), 578–87. doi:10.1043/2010-0613-RAIR.1
- Bergethon, K., Shaw, A. T., Ou, S.-H. I., Katayama, R., Lovly, C. M., McDonald, N. T., ... Iafrate, a J. (2012). ROS1 rearrangements define a unique molecular class of lung cancers. *Journal of Clinical Oncology: Official Journal of the American Society of Clinical Oncology*, 30(8), 863–70. doi:10.1200/JCO.2011.35.6345
- Carvajal, R. D., Antonescu, C. R., Wolchok, J. D., Chapman, P. B., Roman, R.-A., Teitcher, J., ... Schwartz, G. K. (2011). KIT as a therapeutic target in metastatic melanoma. *JAMA: The Journal of the American Medical Association*, 305(22), 2327–34. doi:10.1001/jama.2011.746
- Collisson, E. a., Campbell, J. D., Brooks, A. N., Berger, A. H., Lee, W., Chmielecki, J., ... Sofia, H. J. (2014). Comprehensive molecular profiling of lung adenocarcinoma. *Nature*, 511(7511), 543–50. doi:10.1038/nature13385
- Colombino, M., Capone, M., Lissia, A., Cossu, A., Rubino, C., De Giorgi, V., ... Palmieri, G. (2012). BRAF/NRAS mutation frequencies among primary tumors and metastases in patients with melanoma. *Journal of Clinical Oncology: Official Journal of the American Society of Clinical Oncology*, 30(20), 2522–9. doi:10.1200/JCO.2011.41.2452
- Corless, C. L., Schroeder, A., Griffith, D., Town, A., McGreevey, L., Harrell, P., ... Heinrich, M. C. (2005). PDGFRA mutations in gastrointestinal stromal tumors: frequency, spectrum and in vitro sensitivity to imatinib. *Journal of Clinical Oncology: Official Journal of the American Society of Clinical Oncology*, 23(23), 5357–64. doi:10.1200/JCO.2005.14.068
- Da Rocha Dias, S., Salmonson, T., van Zwieten-Boot, B., Jonsson, B., Marchetti, S., Schellens, J. H. M., ... Pignatti, F. (2013). The European Medicines Agency review of vemurafenib (Zelboraf®) for the treatment of adult patients with BRAF V600 mutation-positive unresectable or metastatic melanoma: Summary of the scientific assessment of the Committee for Medicinal Products for Human. *European Journal of Cancer (Oxford, England: 1990)*, 49(7), 1654–1661. doi:10.1016/j.ejca.2013.01.015
- De Roock, W., Claes, B., Bernasconi, D., De Schutter, J., Biesmans, B., Fountzilias, G., ... Tejpar, S. (2010). Effects of KRAS, BRAF, NRAS, and PIK3CA mutations on

- the efficacy of cetuximab plus chemotherapy in chemotherapy-refractory metastatic colorectal cancer: a retrospective consortium analysis. *The Lancet Oncology*, 11(8), 753–62. doi:10.1016/S1470-2045(10)70130-3
- Douillard, J.-Y., Oliner, K. S., Siena, S., Tabernero, J., Burkes, R., Barugel, M., ... Patterson, S. D. (2013). Panitumumab–FOLFOLX4 Treatment and RAS Mutations in Colorectal Cancer. *New England Journal of Medicine*, 369(11), 1023–1034. doi:10.1056/NEJMoa1305275
- Heinrich, M. C., Corless, C. L., Demetri, G. D., Blanke, C. D., von Mehren, M., Joensuu, H., ... Fletcher, J. a. (2003). Kinase mutations and imatinib response in patients with metastatic gastrointestinal stromal tumor. *Journal of Clinical Oncology : Official Journal of the American Society of Clinical Oncology*, 21(23), 4342–9. doi:10.1200/JCO.2003.04.190
- Heinrich, M. C., Maki, R. G., Corless, C. L., Antonescu, C. R., Harlow, A., Griffith, D., ... Demetri, G. D. (2008). Primary and secondary kinase genotypes correlate with the biological and clinical activity of sunitinib in imatinib-resistant gastrointestinal stromal tumor. *Journal of Clinical Oncology*, 26(33), 5352–5359. doi:10.1200/JCO.2007.15.7461
- Hodi, F. S., Corless, C. L., Giobbie-Hurder, A., Fletcher, J. a., Zhu, M., Marino-Enriquez, A., ... Fisher, D. E. (2013). Imatinib for melanomas harboring mutationally activated or amplified kit arising on mucosal, acral, and chronically sun-damaged skin. *Journal of Clinical Oncology*, 31(26), 3182–3190. doi:10.1200/JCO.2012.47.7836
- Hodis, E., Watson, I. R., Kryukov, G. V, Arold, S. T., Imielinski, M., Theurillat, J.-P., ... Chin, L. (2012). A landscape of driver mutations in melanoma. *Cell*, 150(2), 251–63. doi:10.1016/j.cell.2012.06.024
- Huang, C., Zou, Q., Liu, H., Qiu, B., Li, Q., Lin, Y., Liang, Y. (2020) Management of Non-small cell lung cancer patients with MET exon 14 skipping mutations
- Jakob, J. a., Bassett, R. L., Ng, C. S., Curry, J. L., Joseph, R. W., Alvarado, G. C., ... Davies, M. a. (2012). NRAS mutation status is an independent prognostic factor in metastatic melanoma. *Cancer*, 118(16), 4014–4023. doi:10.1002/cncr.26724
- Kawano, O., Sasaki, H., Endo, K., Suzuki, E., Haneda, H., Yukiue, H., ... Fujii, Y. (2006). PIK3CA mutation status in Japanese lung cancer patients. *Lung Cancer*, 54(2), 209–215. doi:10.1016/j.lungcan.2006.07.006
- Kwak, E. L., Bang, Y.-J., Camidge, D. R., Shaw, A. T., Solomon, B., Maki, R. G., ... Iafrate, A. J. (2010). Anaplastic lymphoma kinase inhibition in non-small-cell lung cancer. *The New England Journal of Medicine*, 363(18), 1693–703. doi:10.1056/NEJMoa1006448
- Lindeman, N. I., Cagle, P. T., Beasley, M. B., Chitale, D. A., Dacic, S., Giaccone, G., ... Ladanyi, M. (2013a). Molecular testing guideline for selection of lung cancer

patients for EGFR and ALK tyrosine kinase inhibitors: guideline from the College of American Pathologists, International Association for the Study of Lung Cancer, and Association for Molecular Patho. *Journal of Thoracic Oncology : Official Publication of the International Association for the Study of Lung Cancer*, 8(7), 823–59. doi:10.1097/JTO.0b013e318290868f

Lindeman, N. I., Cagle, P. T., Beasley, M. B., Chitale, D. A., Dacic, S., Giaccone, G., ... Ladanyi, M. (2013b). Molecular testing guideline for selection of lung cancer patients for EGFR and ALK tyrosine kinase inhibitors: guideline from the College of American Pathologists, International Association for the Study of Lung Cancer, and Association for Molecular Patho. *Archives of Pathology & Laboratory Medicine*, 137(6), 828–60. doi:10.5858/arpa.2012-0720-OA

Lochhead, P., Kuchiba, A., Imamura, Y., Liao, X., Yamauchi, M., Nishihara, R., ... Ogino, S. (2013). Microsatellite instability and BRAF mutation testing in colorectal cancer prognostication. *Journal of the National Cancer Institute*, 105(15), 1151–6. doi:10.1093/jnci/djt173

Miranda, C., Nucifora, M., Molinari, F., Conca, E., Anania, M. C., Bordoni, A., ... Frattini, M. (2012). KRAS and BRAF mutations predict primary resistance to imatinib in gastrointestinal stromal tumors. *Clinical Cancer Research*, 18(6), 1769–1776. doi:10.1158/1078-0432.CCR-11-2230

Mok, T. S., Wu, Y., Thongprasert, S., Yang, C.-H., Chu, D.-T., Saijo, N., ... Fukuoka, M. (2009). Gefitinib or carboplatin-paclitaxel in pulmonary adenocarcinoma. *The New England Journal of Medicine*, 361(10), 947–57. doi:10.1056/NEJMoa0810699

Nazarian, R., Shi, H., Wang, Q., Kong, X., Koya, R. C., Lee, H., ... Lo, R. S. (2010). Melanomas acquire resistance to B-RAF(V600E) inhibition by RTK or N-RAS upregulation. *Nature*, 468(7326), 973–977. doi:10.1038/nature09626

Neumann, J., Horst, D., Kriegl, L., Maatz, S., Engel, J., Jung, A., & Kirchner, T. (2012). A simple immunohistochemical algorithm predicts the risk of distant metastases in right-sided colon cancer. *Histopathology*, 60(3), 416–26. doi:10.1111/j.1365-2559.2011.04126.x

Oxnard, G. R., Miller, V. A., Robson, M. E., Azzoli, C. G., Pao, W., Ladanyi, M., & Arcila, M. E. (2012). Screening for germline EGFR T790M mutations through lung cancer genotyping. *Journal of Thoracic Oncology : Official Publication of the International Association for the Study of Lung Cancer*, 7(6), 1049–52. doi:10.1097/JTO.0b013e318250ed9d

Paik, P. K., Arcila, M. E., Fara, M., Sima, C. S., Miller, V. a., Kris, M. G., ... Riely, G. J. (2011). Clinical characteristics of patients with lung adenocarcinomas harboring BRAF mutations. *Journal of Clinical Oncology*, 29(15), 2046–2051. doi:10.1200/JCO.2010.33.1280

- Pirker, R., Herth, F. J. F., Kerr, K. M., Filipits, M., Taron, M., Gandara, D., ... Stahel, R. (2010). Consensus for EGFR mutation testing in non-small cell lung cancer: results from a European workshop. *Journal of Thoracic Oncology: Official Publication of the International Association for the Study of Lung Cancer*, 5(10), 1706–13. doi:10.1097/JTO.0b013e3181f1c8de
- Rikova, K., Guo, A., Zeng, Q., Possemato, A., Yu, J., Haack, H., ... Comb, M. J. (2007). Global Survey of Phosphotyrosine Signaling Identifies Oncogenic Kinases in Lung Cancer. *Cell*, 131(6), 1190–1203. doi:10.1016/j.cell.2007.11.025
- Shaw, A. T., Kim, D.-W., Nakagawa, K., Seto, T., Crinó, L., Ahn, M.-J., ... Jänne, P. a. (2013). Crizotinib versus chemotherapy in advanced ALK-positive lung cancer. *The New England Journal of Medicine*, 368(25), 2385–94. doi:10.1056/NEJMoa1214886
- Sonobe, M., Manabe, T., Wada, H., & Tanaka, F. (2006). Lung adenocarcinoma harboring mutations in the ERBB2 kinase domain. *The Journal of Molecular Diagnostics : JMD*, 8(3), 351–356. doi:10.2353/jmoldx.2006.050132
- Takeuchi, K., Choi, Y. L., Togashi, Y., Soda, M., Hatano, S., Inamura, K., ... Mano, H. (2009). KIF5B-ALK, a novel fusion oncokinase identified by an immunohistochemistry-based diagnostic system for ALK-positive lung cancer. *Clinical Cancer Research: An Official Journal of the American Association for Cancer Research*, 15(9), 3143–9. doi:10.1158/1078-0432.CCR-08-3248
- TCGA. (2012). Comprehensive genomic characterization of squamous cell lung cancers. *Nature*, 489(7417), 519–25. doi:10.1038/nature11404
- Van den Hurk, K., Balint, B., Toomey, S., O’Leary, P. C., Unwin, L., Sheahan, K., ... Gallagher, W. M. (2015). High-throughput oncogene mutation profiling shows demographic differences in BRAF mutation rates among melanoma patients. *Melanoma Research*. doi:10.1097/CMR.0000000000000149